

Table of contents

\\Projects

FPN

AlessandraPizzuti

sub-07_sess-01

AAHead_Scout_32ch-head-coil
cmrr_mbep2d_bold_pt8_G3_MB2_AP_phy_run1
cmrr_mbep2d_bold_pt8_G3_MB2_AP_phy_run1
localizer_cmrr_mbep2d_bold_1pt8_G2_MB3_AP_run1_ML
mp2rage_iso0.7mm_iPAT3

\\Projects\FPN\AlessandraPizzuti\sub-07_sess-01\AAHead_Scout_32ch-head-coil

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

Properties

| | |
|---|--------------------|
| Prio recon | Off |
| Load images to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load images to stamp segments | Off |
| Load images to graphic segments | On |
| Auto open inline display | Off |
| Auto close inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | Single measurement |

Routine

| | |
|--------------------|--------------------|
| Slab group | 1 |
| Slabs | 1 |
| Dist. factor | 20 % |
| Position | L0.0 A30.0 H0.0 mm |
| Orientation | Sagittal |
| Phase enc. dir. | A >> P |
| Phase oversampling | 0 % |
| Slice oversampling | 0.0 % |
| Slices per slab | 128 |
| FoV read | 260 mm |
| FoV phase | 100.0 % |
| Slice thickness | 1.6 mm |
| TR | 3.25 ms |
| TE | 1.53 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | B1 filter |
| Coil elements | A32 |

Contrast - Common

| | |
|------------|---------|
| TR | 3.25 ms |
| TE | 1.53 ms |
| Flip angle | 16 deg |

Contrast - Dynamic

| | |
|----------------|------------|
| Averages | 1 |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |

Resolution - Common

| | |
|-----------------------|-----------|
| FoV read | 260 mm |
| FoV phase | 100.0 % |
| Slice thickness | 1.6 mm |
| Base resolution | 160 |
| Phase resolution | 100 % |
| Slice resolution | 69 % |
| Phase partial Fourier | 6/8 |
| Slice partial Fourier | 6/8 |
| Trajectory | Cartesian |

Resolution - iPAT

| | |
|------------------|--------|
| PAT mode | GRAPPA |
| Accel. factor PE | 3 |
| Ref. lines PE | 24 |
| Accel. factor 3D | 1 |

Resolution - iPAT

| | |
|---------------------|------------|
| Reference scan mode | Integrated |
|---------------------|------------|

Resolution - Filter Image

| | |
|-------------------|-----|
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | On |
| Unfiltered images | Off |

Resolution - Filter Rawdata

| | |
|-------------------|-----|
| Raw filter | Off |
| Elliptical filter | Off |

Geometry - Common

| | |
|--------------------|--------------------|
| Slab group | 1 |
| Slabs | 1 |
| Dist. factor | 20 % |
| Position | L0.0 A30.0 H0.0 mm |
| Orientation | Sagittal |
| Phase enc. dir. | A >> P |
| Slice oversampling | 0.0 % |
| Slices per slab | 128 |
| FoV read | 260 mm |
| FoV phase | 100.0 % |
| Slice thickness | 1.6 mm |
| TR | 3.25 ms |
| Multi-slice mode | Sequential |
| Series | Ascending |
| Concatenations | 1 |

Geometry - AutoAlign

| | |
|---------------------|--------------------|
| Slab group | 1 |
| Position | L0.0 A30.0 H0.0 mm |
| Orientation | Sagittal |
| Phase enc. dir. | A >> P |
| Initial Position | Isocenter |
| L | 0.0 mm |
| P | 0.0 mm |
| H | 0.0 mm |
| Initial Rotation | 0.00 deg |
| Initial Orientation | Transversal |

Geometry - Tim Planning Suite

| | |
|-------------------|------|
| Set-n-Go Protocol | Off |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System - Miscellaneous

| | |
|---------------------|----------------|
| Positioning mode | REF |
| Table position | H |
| Table position | 0 mm |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| Save uncombined | Off |
| Matrix Optimization | Off |

System - Miscellaneous

| | |
|------------------|---------|
| Coil Select Mode | Default |
|------------------|---------|

System - Adjustments

| | |
|--------------------------|----------|
| B0 Shim mode | Tune up |
| B1 Shim mode | TrueForm |
| Confirm freq. adjustment | Off |
| Assume Dominant Fat | Off |
| Assume Silicone | Off |
| Adjustment Tolerance | Auto |

System - Adjust Volume

| | |
|-------------|-------------|
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| A >> P | 263 mm |
| R >> L | 350 mm |
| F >> H | 350 mm |
| Reset | Off |

System - Tx/Rx

| | |
|---------------------|----------------|
| Frequency 1H | 297.204937 MHz |
| Correction factor | 1 |
| Gain | High |
| Img. Scale Cor. | 1.000 |
| Reset | Off |
| ? Ref. amplitude 1H | 0.000 V |

Physio - PACE

| | |
|----------------|-----|
| Resp. control | Off |
| Concatenations | 1 |

Inline - Common

| | |
|----------------|--------|
| Flip angle | 16 deg |
| Measurements | 1 |
| Time to center | 6.3 s |

Inline - Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Measurements | 1 |
| StdDev | Off |
| Save original images | On |

Inline - MIP

| | |
|----------------------|-----|
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |

Inline - Composing

| | |
|------------------|-----|
| Inline Composing | Off |
| Distortion Corr. | Off |

Inline - MapIt

| | |
|----------------------|---------|
| Save original images | On |
| MapIt | None |
| Flip angle | 16 deg |
| Measurements | 1 |
| Contrasts | 1 |
| TR | 3.25 ms |
| TE | 1.53 ms |

Sequence - Part 1

| | |
|------------------|------------|
| Introduction | On |
| Dimension | 3D |
| Asymmetric echo | Weak |
| Contrasts | 1 |
| Multi-slice mode | Sequential |
| Bandwidth | 540 Hz/Px |

Sequence - Part 2

| | |
|---------------|----------|
| RF pulse type | Fast |
| Gradient mode | Normal |
| Excitation | Non-sel. |
| RF spoiling | On |

Sequence - Nuclei

| | |
|-----------------------|------|
| TX/RX Nucleus | 1H |
| TX/RX delta frequency | 0 Hz |
| TX Nucleus | None |
| TX delta frequency | 0 Hz |
| Coil elements | A32 |

Sequence - Assistant

| | |
|------|-----|
| Mode | Off |
|------|-----|

\\Projects\FPN\AlessandraPizzuti\sub-07_sess-01\cmrr_mbep2d_bold_pt8_G3_MB2_AP_phy_run1

TA: 0:31 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : epfid

Properties

| | |
|---|--------------------|
| Prio recon | Off |
| Load images to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load images to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Auto close inline display | Off |
| Start measurement without further preparation | Off |
| Wait for user to start | Off |
| Start measurements | Single measurement |

Routine

| | |
|--------------------------|--------------------|
| Slice group | 1 |
| Slices | 58 |
| Dist. factor | 0 % |
| Position | L7.1 P4.0 H44.7 mm |
| Orientation | T > C-42.7 > S-6.0 |
| Phase enc. dir. | P >> A |
| AutoAlign | Head > Brain |
| Phase oversampling | 0 % |
| FoV read | 140 mm |
| FoV phase | 97.7 % |
| Slice thickness | 0.80 mm |
| TR | 2000 ms |
| TE | 24.60 ms |
| Multi-band accel. factor | 2 |
| Filter | None |
| Coil elements | A32 |

Contrast - Common

| | |
|-------------------|----------|
| TR | 2000 ms |
| TE | 24.60 ms |
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 69 deg |
| Fat suppr. | Fat sat. |

Contrast - Dynamic

| | |
|-----------------|-------------|
| Averaging mode | Long term |
| Reconstruction | Magn./Phase |
| Measurements | 5 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Resolution - Common

| | |
|-----------------------|---------|
| FoV read | 140 mm |
| FoV phase | 97.7 % |
| Slice thickness | 0.80 mm |
| Base resolution | 176 |
| Phase resolution | 100 % |
| Phase partial Fourier | 6/8 |
| Interpolation | Off |

Resolution - iPAT

| | |
|---------------------|--------|
| PAT mode | GRAPPA |
| Accel. factor PE | 3 |
| Ref. lines PE | 69 |
| Reference scan mode | FLEET |

Resolution - Filter Image

| | |
|-------------------|-----|
| Distortion Corr. | Off |
| Prescan Normalize | Off |

Resolution - Filter Rawdata

| | |
|-------------------|-----|
| Raw filter | Off |
| Elliptical filter | Off |
| Hamming | Off |

Geometry - Common

| | |
|--------------------------|--------------------|
| Slice group | 1 |
| Slices | 58 |
| Dist. factor | 0 % |
| Position | L7.1 P4.0 H44.7 mm |
| Orientation | T > C-42.7 > S-6.0 |
| Phase enc. dir. | P >> A |
| FoV read | 140 mm |
| FoV phase | 97.7 % |
| Slice thickness | 0.80 mm |
| TR | 2000 ms |
| Multi-slice mode | Interleaved |
| Series | Ascending |
| Multi-band accel. factor | 2 |

Geometry - AutoAlign

| | |
|---------------------|--------------------|
| Slice group | 1 |
| Position | L7.1 P4.0 H44.7 mm |
| Orientation | T > C-42.7 > S-6.0 |
| Phase enc. dir. | P >> A |
| AutoAlign | Head > Brain |
| Initial Position | L4.7 P38.9 H23.0 |
| L | 4.7 mm |
| P | 38.9 mm |
| H | 23.0 mm |
| Initial Rotation | 90.20 deg |
| Initial Orientation | C > T |
| C > T | -35.8 |
| > S | -0.8 |

Geometry - Saturation

| | |
|--------------|----------|
| Fat suppr. | Fat sat. |
| Special sat. | None |

Geometry - Tim Planning Suite

| | |
|-------------------|------|
| Set-n-Go Protocol | Off |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System - Miscellaneous

| | |
|---------------------|----------------|
| Positioning mode | FIX |
| Table position | H |
| Table position | 0 mm |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| Matrix Optimization | Off |
| AutoAlign | Head > Brain |
| Coil Select Mode | Default |

System - Adjustments

| | |
|--------------------------|----------|
| B0 Shim mode | Brain |
| B1 Shim mode | TrueForm |
| Confirm freq. adjustment | Off |
| Assume Dominant Fat | Off |
| Assume Silicone | Off |
| Adjustment Tolerance | Auto |

System - Adjust Volume

| | |
|---------------|--------------------|
| ! Position | L5.0 P4.9 H46.3 mm |
| ! Orientation | T > C-42.6 > S-6.9 |
| ! Rotation | -0.32 deg |
| ! A >> P | 135 mm |
| ! R >> L | 149 mm |
| ! F >> H | 42 mm |
| Reset | Off |

System - Tx/Rx

| | |
|---------------------|----------------|
| Frequency 1H | 297.204937 MHz |
| Correction factor | 1 |
| Gain | High |
| Img. Scale Cor. | 1.000 |
| Reset | Off |
| ? Ref. amplitude 1H | 0.000 V |

Physio - Signal1

| | |
|--------------------------|---------|
| 1st Signal/Mode | None |
| TR | 2000 ms |
| Multi-band accel. factor | 2 |

BOLD

| | |
|-------------------------|----------|
| GLM Statistics | Off |
| Dynamic t-maps | Off |
| Ignore meas. at start | 0 |
| Ignore after transition | 0 |
| Model transition states | On |
| Temp. highpass filter | On |
| Threshold | 4.00 |
| Paradigm size | 20 |
| Meas[1] | Baseline |
| Meas[2] | Baseline |
| Meas[3] | Baseline |
| Meas[4] | Baseline |
| Meas[5] | Baseline |
| Meas[6] | Baseline |
| Meas[7] | Baseline |
| Meas[8] | Baseline |
| Meas[9] | Baseline |
| Meas[10] | Baseline |
| Meas[11] | Active |
| Meas[12] | Active |
| Meas[13] | Active |
| Meas[14] | Active |
| Meas[15] | Active |
| Meas[16] | Active |
| Meas[17] | Active |
| Meas[18] | Active |
| Meas[19] | Active |
| Meas[20] | Active |
| Motion correction | Off |
| Spatial filter | Off |
| Measurements | 5 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Sequence - Part 1

| | |
|-------------------|-------------|
| Introduction | On |
| Contrasts | 1 |
| Flow comp. | No |
| Multi-slice mode | Interleaved |
| Free echo spacing | Off |
| Echo spacing | 1.08 ms |
| Bandwidth | 1052 Hz/Px |

Sequence - Part 2

| | |
|---------------|-------|
| EPI factor | 172 |
| Gradient mode | Fast* |
| RF spoiling | Off |

Sequence - Special

| | |
|--------------------------|-----------|
| Excite pulse duration | 4000 us |
| FLEET ref. prep. scans | 0 |
| FLEET ref. min. TR | 0 ms |
| Slice multiplier | 1 |
| EPI noise scans | 0 |
| Single-band images | On |
| MB LeakBlock kernel | Off |
| MB dual kernel | Off |
| MB RF phase scramble | Off |
| SENSE1 coil combine | On |
| Invert RO/PE polarity | On |
| Force Maxwell corr. | Off |
| PF omits higher k-space | Off |
| Disable freq. update | Off |
| Force equal slice timing | Off |
| Online multi-band recon. | Online |
| FFT scale factor | 1.00 |
| Fat saturation FA | 110.0 deg |
| FLEET iPAT ref. FA | 12.0 deg |
| Physio recording | Legacy |
| Triggering scheme | Standard |

\\Projects\FPN\AlessandraPizzuti\sub-07_sess-01\cmrr_mbep2d_bold_pt8_G3_MB2_AP_phy_run1

TA: 10:37 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : epfid

Properties

| | |
|---|--------------------|
| Prio recon | Off |
| Load images to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load images to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Auto close inline display | Off |
| Start measurement without further preparation | Off |
| Wait for user to start | Off |
| Start measurements | Single measurement |

Routine

| | |
|--------------------------|--------------------|
| Slice group | 1 |
| Slices | 58 |
| Dist. factor | 0 % |
| Position | L7.1 P4.0 H44.7 mm |
| Orientation | T > C-42.7 > S-6.0 |
| Phase enc. dir. | P >> A |
| AutoAlign | Head > Brain |
| Phase oversampling | 0 % |
| FoV read | 140 mm |
| FoV phase | 97.7 % |
| Slice thickness | 0.80 mm |
| TR | 2000 ms |
| TE | 24.60 ms |
| Multi-band accel. factor | 2 |
| Filter | None |
| Coil elements | A32 |

Contrast - Common

| | |
|-------------------|----------|
| TR | 2000 ms |
| TE | 24.60 ms |
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 69 deg |
| Fat suppr. | Fat sat. |

Contrast - Dynamic

| | |
|-----------------|-------------|
| Averaging mode | Long term |
| Reconstruction | Magn./Phase |
| Measurements | 308 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Resolution - Common

| | |
|-----------------------|---------|
| FoV read | 140 mm |
| FoV phase | 97.7 % |
| Slice thickness | 0.80 mm |
| Base resolution | 176 |
| Phase resolution | 100 % |
| Phase partial Fourier | 6/8 |
| Interpolation | Off |

Resolution - iPAT

| | |
|---------------------|--------|
| PAT mode | GRAPPA |
| Accel. factor PE | 3 |
| Ref. lines PE | 69 |
| Reference scan mode | FLEET |

Resolution - Filter Image

| | |
|-------------------|-----|
| Distortion Corr. | Off |
| Prescan Normalize | Off |

Resolution - Filter Rawdata

| | |
|-------------------|-----|
| Raw filter | Off |
| Elliptical filter | Off |
| Hamming | Off |

Geometry - Common

| | |
|--------------------------|--------------------|
| Slice group | 1 |
| Slices | 58 |
| Dist. factor | 0 % |
| Position | L7.1 P4.0 H44.7 mm |
| Orientation | T > C-42.7 > S-6.0 |
| Phase enc. dir. | P >> A |
| FoV read | 140 mm |
| FoV phase | 97.7 % |
| Slice thickness | 0.80 mm |
| TR | 2000 ms |
| Multi-slice mode | Interleaved |
| Series | Ascending |
| Multi-band accel. factor | 2 |

Geometry - AutoAlign

| | |
|---------------------|--------------------|
| Slice group | 1 |
| Position | L7.1 P4.0 H44.7 mm |
| Orientation | T > C-42.7 > S-6.0 |
| Phase enc. dir. | P >> A |
| AutoAlign | Head > Brain |
| Initial Position | L4.7 P38.9 H23.0 |
| L | 4.7 mm |
| P | 38.9 mm |
| H | 23.0 mm |
| Initial Rotation | 90.20 deg |
| Initial Orientation | C > T |
| C > T | -35.8 |
| > S | -0.8 |

Geometry - Saturation

| | |
|--------------|----------|
| Fat suppr. | Fat sat. |
| Special sat. | None |

Geometry - Tim Planning Suite

| | |
|-------------------|------|
| Set-n-Go Protocol | Off |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System - Miscellaneous

| | |
|---------------------|----------------|
| Positioning mode | FIX |
| Table position | H |
| Table position | 0 mm |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| Matrix Optimization | Off |
| AutoAlign | Head > Brain |
| Coil Select Mode | Default |

System - Adjustments

| | |
|--------------------------|----------|
| B0 Shim mode | Brain |
| B1 Shim mode | TrueForm |
| Confirm freq. adjustment | Off |
| Assume Dominant Fat | Off |
| Assume Silicone | Off |
| Adjustment Tolerance | Auto |

System - Adjust Volume

| | |
|---------------|--------------------|
| ! Position | L5.0 P4.9 H46.3 mm |
| ! Orientation | T > C-42.6 > S-6.9 |
| ! Rotation | -0.32 deg |
| ! A >> P | 135 mm |
| ! R >> L | 149 mm |
| ! F >> H | 42 mm |
| Reset | Off |

System - Tx/Rx

| | |
|---------------------|----------------|
| Frequency 1H | 297.204937 MHz |
| Correction factor | 1 |
| Gain | High |
| Img. Scale Cor. | 1.000 |
| Reset | Off |
| ? Ref. amplitude 1H | 0.000 V |

Physio - Signal1

| | |
|--------------------------|---------|
| 1st Signal/Mode | None |
| TR | 2000 ms |
| Multi-band accel. factor | 2 |

BOLD

| | |
|-------------------------|----------|
| GLM Statistics | Off |
| Dynamic t-maps | Off |
| Ignore meas. at start | 0 |
| Ignore after transition | 0 |
| Model transition states | On |
| Temp. highpass filter | On |
| Threshold | 4.00 |
| Paradigm size | 20 |
| Meas[1] | Baseline |
| Meas[2] | Baseline |
| Meas[3] | Baseline |
| Meas[4] | Baseline |
| Meas[5] | Baseline |
| Meas[6] | Baseline |
| Meas[7] | Baseline |
| Meas[8] | Baseline |
| Meas[9] | Baseline |
| Meas[10] | Baseline |
| Meas[11] | Active |
| Meas[12] | Active |
| Meas[13] | Active |
| Meas[14] | Active |
| Meas[15] | Active |
| Meas[16] | Active |
| Meas[17] | Active |
| Meas[18] | Active |
| Meas[19] | Active |
| Meas[20] | Active |
| Motion correction | Off |
| Spatial filter | Off |
| Measurements | 308 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Sequence - Part 1

| | |
|-------------------|-------------|
| Introduction | On |
| Contrasts | 1 |
| Flow comp. | No |
| Multi-slice mode | Interleaved |
| Free echo spacing | Off |
| Echo spacing | 1.08 ms |
| Bandwidth | 1052 Hz/Px |

Sequence - Part 2

| | |
|---------------|-------|
| EPI factor | 172 |
| Gradient mode | Fast* |
| RF spoiling | Off |

Sequence - Special

| | |
|--------------------------|-----------|
| Excite pulse duration | 4000 us |
| FLEET ref. prep. scans | 0 |
| FLEET ref. min. TR | 0 ms |
| Slice multiplier | 1 |
| EPI noise scans | 0 |
| Single-band images | On |
| MB LeakBlock kernel | Off |
| MB dual kernel | Off |
| MB RF phase scramble | Off |
| SENSE1 coil combine | On |
| Invert RO/PE polarity | Off |
| Force Maxwell corr. | Off |
| PF omits higher k-space | Off |
| Disable freq. update | Off |
| Force equal slice timing | Off |
| Online multi-band recon. | Online |
| FFT scale factor | 1.00 |
| Fat saturation FA | 110.0 deg |
| FLEET iPAT ref. FA | 12.0 deg |
| Physio recording | Legacy |
| Triggering scheme | Standard |

\\Projects\FPN\AlessandraPizzuti\sub-07_sess-01\localizer_cmrr_mbep2d_bold_1pt8_G2_MB3_AP_run1_ML

TA: 4:48 PM: FIX Voxel size: 1.8×1.8×1.8 mmPAT: 2 Rel. SNR: 1.00 : epfid

Properties

| | |
|---|--------------------|
| Prio recon | Off |
| Load images to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load images to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Auto close inline display | Off |
| Start measurement without further preparation | Off |
| Wait for user to start | Off |
| Start measurements | Single measurement |

Routine

| | |
|--------------------------|---------------------|
| Slice group | 1 |
| Slices | 48 |
| Dist. factor | 0 % |
| Position | R0.6 A20.9 H41.1 mm |
| Orientation | T > C-23.9 > S-6.0 |
| Phase enc. dir. | A >> P |
| AutoAlign | --- |
| Phase oversampling | 0 % |
| FoV read | 180 mm |
| FoV phase | 100.0 % |
| Slice thickness | 1.80 mm |
| TR | 1000 ms |
| TE | 23.00 ms |
| Multi-band accel. factor | 3 |
| Filter | None |
| Coil elements | A32 |

Contrast - Common

| | |
|-------------------|----------|
| TR | 1000 ms |
| TE | 23.00 ms |
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 54 deg |
| Fat suppr. | Fat sat. |

Contrast - Dynamic

| | |
|-----------------|-----------|
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 270 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Resolution - Common

| | |
|-----------------------|---------|
| FoV read | 180 mm |
| FoV phase | 100.0 % |
| Slice thickness | 1.80 mm |
| Base resolution | 100 |
| Phase resolution | 100 % |
| Phase partial Fourier | 6/8 |
| Interpolation | Off |

Resolution - iPAT

| | |
|------------------|--------|
| PAT mode | GRAPPA |
| Accel. factor PE | 2 |
| Ref. lines PE | 24 |

Resolution - iPAT

| | |
|---------------------|-------------|
| Reference scan mode | Single-shot |
|---------------------|-------------|

Resolution - Filter Image

| | |
|-------------------|-----|
| Distortion Corr. | Off |
| Prescan Normalize | Off |

Resolution - Filter Rawdata

| | |
|-------------------|-----|
| Raw filter | Off |
| Elliptical filter | Off |
| Hamming | Off |

Geometry - Common

| | |
|--------------------------|---------------------|
| Slice group | 1 |
| Slices | 48 |
| Dist. factor | 0 % |
| Position | R0.6 A20.9 H41.1 mm |
| Orientation | T > C-23.9 > S-6.0 |
| Phase enc. dir. | A >> P |
| FoV read | 180 mm |
| FoV phase | 100.0 % |
| Slice thickness | 1.80 mm |
| TR | 1000 ms |
| Multi-slice mode | Interleaved |
| Series | Ascending |
| Multi-band accel. factor | 3 |

Geometry - AutoAlign

| | |
|---------------------|---------------------|
| Slice group | 1 |
| Position | R0.6 A20.9 H41.1 mm |
| Orientation | T > C-23.9 > S-6.0 |
| Phase enc. dir. | A >> P |
| AutoAlign | --- |
| Initial Position | R0.6 A20.9 H41.1 |
| R | 0.6 mm |
| A | 20.9 mm |
| H | 41.1 mm |
| Initial Rotation | 6.64 deg |
| Initial Orientation | T > C |
| T > C | -23.9 |
| > S | -6.0 |

Geometry - Saturation

| | |
|--------------|----------|
| Fat suppr. | Fat sat. |
| Special sat. | None |

Geometry - Tim Planning Suite

| | |
|-------------------|------|
| Set-n-Go Protocol | Off |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System - Miscellaneous

| | |
|------------------|-----------|
| Positioning mode | FIX |
| Table position | H |
| Table position | 0 mm |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |

System - Miscellaneous

| | |
|---------------------|----------------|
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| Matrix Optimization | Off |
| AutoAlign | --- |
| Coil Select Mode | Default |

System - Adjustments

| | |
|--------------------------|----------|
| B0 Shim mode | Tune up |
| B1 Shim mode | TrueForm |
| Confirm freq. adjustment | Off |
| Assume Dominant Fat | Off |
| Assume Silicone | Off |
| Adjustment Tolerance | Auto |

System - Adjust Volume

| | |
|---------------|---------------------|
| ! Position | R2.3 A21.4 H18.6 mm |
| ! Orientation | T > C-6.6 |
| ! Rotation | 0.27 deg |
| ! A >> P | 191 mm |
| ! R >> L | 176 mm |
| ! F >> H | 96 mm |
| Reset | Off |

System - Tx/Rx

| | |
|---------------------|----------------|
| Frequency 1H | 297.204937 MHz |
| Correction factor | 1 |
| Gain | High |
| Img. Scale Cor. | 1.000 |
| Reset | Off |
| ? Ref. amplitude 1H | 0.000 V |

Physio - Signal1

| | |
|--------------------------|---------|
| 1st Signal/Mode | None |
| TR | 1000 ms |
| Multi-band accel. factor | 3 |

BOLD

| | |
|-------------------------|----------|
| GLM Statistics | On |
| Dynamic t-maps | On |
| Ignore meas. at start | 0 |
| Ignore after transition | 0 |
| Model transition states | On |
| Temp. highpass filter | On |
| Threshold | 4.00 |
| Paradigm size | 20 |
| Meas[1] | Baseline |
| Meas[2] | Baseline |
| Meas[3] | Baseline |
| Meas[4] | Baseline |
| Meas[5] | Baseline |
| Meas[6] | Baseline |
| Meas[7] | Baseline |
| Meas[8] | Baseline |
| Meas[9] | Baseline |
| Meas[10] | Baseline |
| Meas[11] | Active |
| Meas[12] | Active |
| Meas[13] | Active |
| Meas[14] | Active |
| Meas[15] | Active |
| Meas[16] | Active |
| Meas[17] | Active |
| Meas[18] | Active |
| Meas[19] | Active |

BOLD

| | |
|-------------------|--------|
| Meas[20] | Active |
| Motion correction | Off |
| Spatial filter | Off |
| Measurements | 270 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Sequence - Part 1

| | |
|-------------------|-------------|
| Introduction | On |
| Contrasts | 1 |
| Flow comp. | No |
| Multi-slice mode | Interleaved |
| Free echo spacing | Off |
| Echo spacing | 1.02 ms |
| Bandwidth | 1086 Hz/Px |

Sequence - Part 2

| | |
|---------------|------|
| EPI factor | 100 |
| Gradient mode | Fast |
| RF spoiling | Off |

Sequence - Special

| | |
|--------------------------|-----------|
| Excite pulse duration | 4000 us |
| Slice multiplier | 1 |
| EPI noise scans | 0 |
| Single-band images | On |
| MB LeakBlock kernel | Off |
| MB dual kernel | Off |
| MB RF phase scramble | On |
| SENSE1 coil combine | On |
| Invert RO/PE polarity | Off |
| PF omits higher k-space | Off |
| Disable freq. update | Off |
| Force equal slice timing | Off |
| Online multi-band recon. | Online |
| FFT scale factor | 0.20 |
| Fat saturation FA | 110.0 deg |
| Physio recording | Legacy |
| Triggering scheme | Standard |

\\Projects\FPN\AlessandraPizzuti\sub-07_sess-01\mp2rage_iso0.7mm_iPAT3

TA: 8:02 PM: REF Voxel size: 0.7×0.7×0.7 mmPAT: 3 Rel. SNR: 1.00 : tfl

Properties

| | |
|---|--------------------|
| Prio recon | Off |
| Load images to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load images to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Auto close inline display | Off |
| Start measurement without further preparation | Off |
| Wait for user to start | Off |
| Start measurements | Single measurement |

Routine

| | |
|--------------------|---------------------|
| Slab group | 1 |
| Slabs | 1 |
| Dist. factor | 50 % |
| Position | R0.9 A31.6 F19.7 mm |
| Orientation | S > C-0.2 |
| Phase enc. dir. | A >> P |
| AutoAlign | --- |
| Phase oversampling | 0 % |
| Slice oversampling | 0.0 % |
| Slices per slab | 240 |
| FoV read | 224 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.70 mm |
| TR | 5000.0 ms |
| TE | 2.47 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | A32 |

Contrast - Common

| | |
|-------------------|---------------------|
| TR | 5000.0 ms |
| TE | 2.47 ms |
| Magn. preparation | Non-sel. IR |
| TI 1 | 900 ms |
| TI 2 | 2750 ms |
| Flip angle 1 | 5.0 deg |
| Flip angle 2 | 3.0 deg |
| Fat suppr. | Water excit. normal |
| Water suppr. | None |

Contrast - Dynamic

| | |
|-----------------|------------------|
| Averages | 1 |
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution - Common

| | |
|-----------------------|---------|
| FoV read | 224 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.70 mm |
| Base resolution | 320 |
| Phase resolution | 100 % |
| Slice resolution | 100 % |
| Phase partial Fourier | 6/8 |

Resolution - Common

| | |
|-----------------------|-----|
| Slice partial Fourier | Off |
| Interpolation | Off |

Resolution - iPAT

| | |
|---------------------|------------|
| PAT mode | GRAPPA |
| Accel. factor PE | 3 |
| Ref. lines PE | 24 |
| Accel. factor 3D | 1 |
| Reference scan mode | Integrated |

Resolution - Filter Image

| | |
|-------------------|-----|
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |

Resolution - Filter Rawdata

| | |
|-------------------|-----|
| Raw filter | Off |
| Elliptical filter | Off |

Geometry - Common

| | |
|--------------------|---------------------|
| Slab group | 1 |
| Slabs | 1 |
| Dist. factor | 50 % |
| Position | R0.9 A31.6 F19.7 mm |
| Orientation | S > C-0.2 |
| Phase enc. dir. | A >> P |
| Slice oversampling | 0.0 % |
| Slices per slab | 240 |
| FoV read | 224 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.70 mm |
| TR | 5000.0 ms |
| Multi-slice mode | Single shot |
| Series | Interleaved |
| Concatenations | 1 |

Geometry - AutoAlign

| | |
|---------------------|---------------------|
| Slab group | 1 |
| Position | R0.9 A31.6 F19.7 mm |
| Orientation | S > C-0.2 |
| Phase enc. dir. | A >> P |
| AutoAlign | --- |
| Initial Position | R0.9 A31.6 F19.7 |
| R | 0.9 mm |
| A | 31.6 mm |
| F | 19.7 mm |
| Initial Rotation | 6.19 deg |
| Initial Orientation | S > C |
| S > C | -0.2 |
| > T | 0.0 |

Geometry - Navigator**Geometry - Tim Planning Suite**

| | |
|-------------------|------|
| Set-n-Go Protocol | Off |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System - Miscellaneous

| | |
|---------------------|------------------|
| Positioning mode | REF |
| Table position | H |
| Table position | 0 mm |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Adaptive Combine |
| Save uncombined | Off |
| Matrix Optimization | Off |
| AutoAlign | --- |
| Coil Select Mode | Off - All |

System - Adjustments

| | |
|--------------------------|----------|
| B0 Shim mode | Tune up |
| B1 Shim mode | TrueForm |
| Confirm freq. adjustment | Off |
| Assume Dominant Fat | Off |
| Assume Silicone | Off |
| Adjustment Tolerance | Auto |

System - Adjust Volume

| | |
|---------------|---------------------|
| ! Position | R2.3 A21.4 H18.6 mm |
| ! Orientation | T > C-6.6 |
| ! Rotation | 0.27 deg |
| ! A >> P | 191 mm |
| ! R >> L | 176 mm |
| ! F >> H | 96 mm |
| Reset | Off |

System - Tx/Rx

| | |
|---------------------|----------------|
| Frequency 1H | 297.204937 MHz |
| Correction factor | 1 |
| Gain | High |
| Img. Scale Cor. | 1.000 |
| Reset | Off |
| ? Ref. amplitude 1H | 0.000 V |

Physio - Signal1

| | |
|-----------------|-----------|
| 1st Signal/Mode | None |
| TR | 5000.0 ms |
| Concatenations | 1 |

Physio - Cardiac

| | |
|-------------------|---------------------|
| Magn. preparation | Non-sel. IR |
| TI 1 | 900 ms |
| TI 2 | 2750 ms |
| Fat suppr. | Water excit. normal |
| Dark blood | Off |
| FoV read | 224 mm |
| FoV phase | 100.0 % |
| Phase resolution | 100 % |

Physio - PACE

| | |
|----------------|-----|
| Resp. control | Off |
| Concatenations | 1 |

Inline - Common

| | |
|----------------------|-----|
| Subtract | Off |
| Measurements | 1 |
| StdDev | Off |
| Save original images | On |

Inline - MIP

| | |
|----------------------|-----|
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |

Inline - Composing

| | |
|------------------|-----|
| Inline Composing | Off |
| Distortion Corr. | Off |

Inline - MapIt

| | |
|----------------------|-----------|
| Save original images | On |
| MapIt | T1 map |
| Flip angle 1 | 5.0 deg |
| Flip angle 2 | 3.0 deg |
| Measurements | 1 |
| TR | 5000.0 ms |
| TE | 2.47 ms |

Sequence - Part 1

| | |
|---------------------|-------------|
| Introduction | On |
| Dimension | 3D |
| Elliptical scanning | Off |
| Reordering | Linear |
| Asymmetric echo | Allowed |
| Flow comp. | No |
| Multi-slice mode | Single shot |
| Echo spacing | 6.9 ms |
| Bandwidth | 250 Hz/Px |

Sequence - Part 2

| | |
|-------------------------|----------|
| RF pulse type | Fast |
| Gradient mode | Fast |
| Excitation | Non-sel. |
| RF spoiling | On |
| Incr. Gradient spoiling | Off |
| Turbo factor | 240 |

Sequence - Nuclei

| | |
|-----------------------|------|
| TX/RX Nucleus | 1H |
| TX/RX delta frequency | 0 Hz |
| TX Nucleus | None |
| TX delta frequency | 0 Hz |
| Coil elements | A32 |

Sequence - Assistant

| | |
|------|-----|
| Mode | Off |
|------|-----|